BIOLUX Advance

Art. No. 52-03000

-  **Bedienungsanleitung**
-  **Operating Instructions**
-  **Mode d'emploi**
-  **Instrucciones de uso**

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General Information

About this Instruction Manual

These operating instructions are to be considered a component of the device. Please read the safety instructions and the operating instructions carefully before use. Keep these instructions for renewed use at a later date. When the device is sold or given to someone else, the instruction manual must be provided to the new owner/user of the product.

DANGER!

You will find this symbol in front of every section of text which deals with the risk of severe injury or even death in the event of improper use.

CAREFUL!

You will find this symbol before every section of text that deals with the risk of minor to severe injuries resulting from improper use.

NOTE!

You will find this symbol in front of every section of text which deals with the risk of damage to property or the environment.

Intended Use

This product is intended only for private use. It was developed for the magnified display of things in nature.

General Warnings

The device is intended only for indoor use.

DANGER!

Tools with sharp edges and points are often used when working with this device. Because there is a risk of injury from such tools, store this device and all tools and accessories in a location that is out of the reach of children.

This product contains small parts which could be swallowed by children. There is a RISK OF CHOKING.

This device contains electronic components which operate via a power source (batteries). Children should only use the device under adult supervision. Only use the device as described in the manual; otherwise, you run the risk of an electric shock.

Do not expose the device to high temperatures. Use only the recommended batteries. Do not short-circuit the device or batteries, or throw them into a fire. Excessive heat or improper handling could trigger a short-circuit, a fire or an explosion.

Leaking battery acid can lead to chemical burns. Avoid contact of battery acid with skin, eyes and mucous membranes. In the event of contact, rinse the affected region immediately with a plenty of water and seek medical attention.

Never bend, pinch or pull the power and connecting cables, extensions and adapters. Protect the cables from sharp edges and heat.

Before operating, check the device, cables and connections for damage. Never use a damaged unit or a unit with damaged power cables. Damaged parts must be exchanged immediately by an authorised service centre.

Children should only use the device under adult supervision. Keep packaging material, like plastic bags and rubber bands, out of the reach of children, as these materials pose a choking hazard.

CAREFUL!

Children must not have access to the included chemicals and liquids. Do not drink the chemicals. Wash hands thoroughly with running water after use. In the event of contact with the eyes or mouth, rinse thoroughly with water. In the event of pain, contact a doctor immediately and take the substances with you.

Before cleaning the device, disconnect it from the power supply by removing the plug or batteries.

! NOTE!

Use only the recommended batteries. Always replace weak or empty batteries with a new, complete set of batteries at full capacity. Do not use batteries from different brands or with different capacities. The batteries should be removed from the unit if it has not been used for a long time.

Never recharge normal, non-rechargeable batteries. This could lead to explosion during the charging process.

Do not disassemble the device. In the event of a defect, please contact your dealer. The dealer will contact the Service Centre and can send the device in to be repaired, if necessary.

Operating instructions

Parts overview (Fig. 1+2):

- ❶ 5x WF Eyepiece
- ❷ USB Interface
- ❸ Barlow lens
- ❹ Built-in digital camera
- ❺ Eyepiece holder
- ❻ Microscope head
- ❼ Objective-revolver
- ❽ Microscope stage
- ❾ Focus wheel
- ❿ LED lighting (transmitted light)
- ⓫ Preparation clamps
- ⓬ Microscope base
- ⓭ Photomizer software
- ⓮ Battery compartment
- ⓯ 5 Slides, 10 Covering glasses and 5 preparations in a plastic box
- ⓰ Dimmer (transmitted light)
- ⓱ Dimmer (direct light)
- ⓲ Accessory tray
- ⓳ Colour filter disc
- ⓴ LED lighting (direct light)
- ⓵ Preparation needles
- ⓶ Microscope tools:
 - a) Pipette; b) Tweezers
- ⓷ Hatchery
- ⓸ MicroCut
- ⓹ Specimens:
 - a) Yeast; b) Gum media (specimen

inclusion medium); c) Sea salt;

d) Brine shrimp eggs

- ⓺ USB connection cable

1. General/Location

First you must make sure that your microscope is on a stable and solid surface.

2. Electric LED lighting with dimmer

The microscope has two independently adjustable lighting units for upper and lower lighting that can be controlled using the dimmer wheel (Fig. 1, 16+17). The transmitted light unit (Fig. 1, 16) is used for transparent specimens on holders/slides. To view solid non-transparent specimens use the direct light unit (Fig. 1, 20). Simultaneous use of both direct and transmitted light is only reasonable if the specimen is semi-transparent. This operating mode is not recommended for transmitted light specimens on slides as it may cause reflection on the slide.

The recommended batteries must be inserted into the battery compartment (Fig. 1, 14) located in the base plate before use.

Remove the battery compartment cover by pressing lightly on the cover clip. Insert the batteries into the holder.

**NOTE!**

Pay attention to the polarity as marked on the battery holder terminals.

As your device has continuously controllable lighting (dimmer) optimal illumination of the object to be viewed is guaranteed.

3. Colour filter disc

The colour filter under the microscope table (Fig. 1, 19) aids in viewing very bright and transparent objects. Just select the right colour for the specimen in question. The components of colourless/transparent objects (e.g. starch particles, single-cell-specimens) can thus be better recognised.

4. Microscope setup

Begin every observation with the lowest magnification.

Place the microscope's table (Fig. 1, 8) with the focus knob (Fig 1, 9) into the lowest position and rotate the objective revolver (Fig. 1, 7) until it locks on the lowest magnification (4x).

**NOTE:**

First, drive the Microscope's table (Fig. 1, 8) to its lowest position before changing the objective in order to prevent damage.

Insert the 5x eyepiece (Fig. 4, 1) in the Barlow lens (Fig. 4, 3).

Take care, that the Barlow lens is inserted completely in the monocular head (Fig. 4, 5).

5. Observation

After you have set up the microscope with the corresponding illumination, the following principles are important:

Begin each observation with a simple observation at lowest magnification, so that the centre and position of the object to be viewed is in focus.

The higher the magnification the more light is required for good picture quality.

Place a slide with permanent with permanent specimen (Fig. 3, 15) right under the objective lens on the microscope table (Fig. 3, 8) and secure it in place with the clamps (Fig. 3, 11). The specimen must be precisely sited in the centre of the lighting (Fig. 3, 10). If it

is not, move the slide to the left and right as well as back and forth until it is centered over the lighting.

Look through the eyepiece (Fig. 1, 1) and turn carefully the focusing wheel (Fig. 1, 9) until you can see a sharp picture.

Now you can get a higher magnification, while you pull out slowly the Barlow lens (Fig. 5, 3) of the monocular barrel (Fig. 5, 5). With the nearly entirely pulled out Barlow lens the magnification is raised to 2x.

**TIP:**

Depending on the preparation higher magnifications do not always lead to better pictures.

With changing magnification (objective lens changes, pulling out of the Barlow lens) the sharpness of the image must be newly defined by turning the focusing wheel (Fig. 1, 9).

**NOTE:**

Please be very careful when doing this. When you move the mechanical plate upwards too fast the objective lens and the slide can touch and become damaged.

6. Viewed Object – condition and preparation


6.1. Condition

With the Barlow lens nearly fully extended magnification can be doubled. Both transparent and non-transparent specimens can be examined with this microscope, which is a direct as well as transmitted light model. If opaque specimens are examined - such as small animals, plant parts, tissue, stone and so on - the light is reflected from the specimen through the lens and eyepiece, where it is magnified, to the eye (reflected light principle, switch position I). If opaque specimens are examined the light from below goes through the specimen, lens and eyepiece to the eye and is magnified en route (direct light principle, switch position II). Many small water organisms, plant parts and the finest animal components have natural transparent characteristics others need to be prepared accordingly. Preparation can take many forms such as the soaking in appropriate substances (medium) to make the object opaque. Objects can also be finely sliced (hand cut, micro-cut). The following sections explain these procedures in detail.

6.2. Creation of thin preparation cuts

Specimens should be sliced as thin as possible, as stated before. A little wax or

paraffin is needed to achieve the best results. A candle can be used for the purpose. The wax is put in a bowl and heated over a flame. The specimen is then dipped several times in the liquid wax. The wax is finally allowed to harden. Use a MicroCut (Fig 2, 24) or knife/scalpel (carefully) to make very thin slices of the object in its wax casing.

These slices are then laid on a glass slide and covered with another.

6.3. Creation of an own preparation

Put the object to be observed on a glass slide and add a drop of distilled water using a pipette (Fig. 8, 22a) onto the object (Fig. 8).

Set a cover glass (available in most well stocked hobby shops) perpendicularly at the edge of the water drop, so that the water runs along the cover glass edge (Fig. 8). Now lower the cover glass slowly over the water drop.

TIP:

The gum medium supplied (Fig 2, 25b) is used to make permanent slide cultures. Add it instead of distilled water. The gum medium hardens so that the specimen is permanently affixed to its slide.

7. Experiments

If you have made yourself familiar with the microscope already, you can accomplish the following experiments and observe the results under your microscope.

7.1. Newspaper print

Objects:

1. A small piece of paper from a newspaper with parts of a picture and some letters
2. A similar piece of paper from an illustrated magazine

Use your microscope at the lowest magnification and use the preparation of the daily paper. The letters seen are broken up, because the newspaper is printed on raw, inferior paper. Letters of the magazines appear smoother and more complete. The picture of the daily paper consists of many small points, which appear somewhat dirty. The pixels (raster points) of the magazine appear sharply.

7.2. Textile fibres

Objects and accessories:

1. Threads of different textiles: Cotton, line, wool, silk, Celanese, nylon etc..
2. Two needles

Each thread is put on a glass slide and frayed with the help of the two needles. The threads are dampened and covered with a cover glass. The microscope is adjusted to a low magnification. Cotton staples are of vegetable origin and look under the microscope like a flat, turned volume. The fibres are thicker and rounder at the edges than in the centre. Cotton staples consist primary of long, collapsed tubes. Linen fibres are also vegetable origin; they are round and run in straight lines direction. The fibres shine like silk and exhibit countless swellings at the fibre pipe. Silk is animal origin and consists of solid fibres of smaller diameter contrary to the hollow vegetable fibres. Each fibre is smooth and even moderate and has the appearance of a small glass rod. Wool fibres are also animal origin; the surface consists of overlapping cases, which appear broken and wavy. If it is possible, compare wool fibres of different weaving mills. Consider thereby the different appearance of the fibres. Experts can determine from it the country of origin of wool. Celanese is like the name says, artificially manufactured by a long chemical process. All fibres show hard, dark lines on

the smooth, shining surface. The fibres ripple themselves/crinkle after drying in the same condition. Observe the thing in common and differences.

7.3. Brine shrimps

Accessories:

1. Brine shrimp eggs (Fig 2, 25d)
2. Sea salt (Fig 2, 25c)
3. Hatchery (Fig 2, 23)
4. Yeast (Fig 2, 25a)

CAREFUL!

Eggs and brine shrimps are not fit for human consumption.

9.3.1 The lifecycle of a brine shrimp

The brine shrimp or *artemia salina* to scientists has an unusual and interesting lifecycle. The female's eggs are hatched without any male shrimp having to fertilise them. The resultant baby shrimps are all female. Under unusual circumstances such as when a swamp is drained the eggs may produce male shrimps. These males fertilise the female's eggs, resulting in a specific type of eggs. These are called winter eggs and have a thick shell as protection. They're pretty rugged and can survive the swamp or lake drying out causing the death of the entire shrimp population for up to a decade in a form of hibernation. The

eggs hatch once again as soon as the correct ambient conditions have been obtained. The eggs supplied (Fig 2, 25d) are of this type.

7.3.2. Hatching of the brine shrimp

To hatch the shrimp it is essential to first have a saline solution suited to the shrimp's needs. Fill half a litre of rain- or fresh-water in a container. Let it stand for about thirty hours. As water evaporates over time it's a good idea to have a second container of such water left standing for thirty-six hours. Once stood for this length of time pour half of the sea salt supplied into one of the containers (Fig 2, 25c) and stir until it has dissolved. Then pour some of it into the hatchery (Fig 2, 23). Add a few eggs and close the lid. Put it somewhere with plenty of light but not in the direct sun. The temperature should be approximately 25° C. The shrimps will hatch in two or three days at this temperature. Should any water evaporate during this time replace it from the second container.

7.3.3. The brine shrimp under the microscope

What comes out of the egg is known as a nauplius larva. Use the pipette (Fig 2, 22a) to put some of them on a slide for examination. They will move in the solution using their hair like limbs. Remove a few daily from the container for examination under the

microscope. If you do so and save the pictures made with the MicrOcular you will then have a seamless record of the shrimp's lifecycle. You can remove the upper lid of the hatchery and put the whole thing under the microscope. The larvae will mature in six to ten weeks depending on ambient temperature. You will soon have bred an entire generation of brine shrimps that constantly reproduce.

7.3.4. Feeding your brine shrimps

To keep them alive brine shrimps must be fed occasionally. This must be done carefully as overfeeding causes the water to stagnate and poison the shrimps. Feeding is best done with dry powdered yeast (Fig 2, 25a). Give them a little every other day. If the water darkens this signifies it is stagnating. If so remove the shrimps and put them in a fresh saline solution.

8. Installation and use of the software

8.1. Software and installation information

A software CD (Fig. 2, 13) is included with your microscope. The software and drivers on this CD must be installed on your computer in order to use the integrated camera. Once you connect the microscope over the USB cable (Fig. 2, 26) to your computer after installation, you can view pictures on and save them to your computer using Photomizer. To install the software and driver correctly, simply follow the installation steps.

8.2. Software installation

1. **Important:** Before inserting the CD, first plug the USB cable (Fig. 2, 26) into the USB port on your PC. Windows will now recognize the new device, and will indicate this in a notification window. Now please click on "Cancel".
2. Now insert the CD-ROM that came in your package into the CD/DVD drive of your computer. The installation menu starts automatically. If it does not, go to the Windows explorer and select the CD/DVD drive (most of the time, it's the "D" drive, but it can have another letter). From there, start the file "autorun.exe" by double-clicking with the left mouse button.

8.2.1. Driver installation

To install the driver software, click on the menu point "Install Driver" with the mouse cursor. Then follow the installation program instructions. During the software installation the correct driver for your operating system will be installed automatically. No manual input is needed. In rare cases the device may not be recognised by your computer. As a rule you need then only install the driver manually from the CD. If this fails please refer to the troubleshooting chapter that follows.

8.2.2. Installation of the Image Editing Software Photomizer

The image editing software "Photomizer" is located on the software CD. You can edit your pictures here.

1. To install the software, click once on the menu point "Install PHOTOMIZER" with the left mouse button.
2. The Photomizer Software requires Microsoft .NET Framework 4.0, which it will install if it is not already on your system. If it is already on your system, you may skip to step 5.
3. In the welcome window, please accept the Microsoft license agreement, and then click "Install". The installation can take a few minutes.
4. As soon as everything is installed, click "Finish".

- Now you will be presented with a choice, in which you can choose your language. Make your selection and confirm it by clicking „OK“.
- When you see “Welcome”, click on „Next“.
- In the next window, you will be asked for the „Destination Folder“. Here, just click on „Next“.
- Now, the window with the Setup Status will appear – here, a progress bar will inform you about the ongoing installation. This process can take a few minutes.
- The window “Photomizer is being installed” appears. Click on “Finish”. The installation ends.

**TIP:**

To use the microscope camera long-term we recommend it always be connected to the same USB port.

9. Using the integrated camera

9.1 Preparation

- Slide a specimen under your microscope and focus on it.
- Start your PC if you haven't yet and connect The microscope to the USB port of your computer.

9.2 Showing and saving images taken with the microscope camera on your PC

- Start the Photomizer Software.
- Click „Open camera“
- In the event that you have connected more than one unit, you can choose the desired unit in the subsequent selection. Click on "USB 2.0 Webcam". If only one device is connected, this step is omitted.
- The camera image should now be visible on your screen. Focus the microscope image.
- Click „Capture“ to record an image. It will then be shown on the right.
- Click on it to select it and then click „Transfer image“.
- This takes you to the Photomizer software.
- File - Save as

9.3 The Photomizer Software

If you need help in the „Photomizer“ programme click „?“ and then „Open help“. If you need further assistance please visit the maker's home page at www.photomizer.net/bresser

Care and maintenance

Separate the device from the power supply before cleaning (remove batteries)!

Only use a dry cloth to clean the exterior of the device.

**NOTE:**

Do not use any cleaning fluid to avoid damaging the electronics.

Clean the lenses (eyepiece and/or lens) with a soft and lint-free cloth only (e.g. microfibre).

**NOTE:**

Do not apply excess pressure to the cloth so as to avoid scratching the lenses.

To remove more stubborn dirt, moisten the cleaning cloth with an eyeglass-cleaning solution and wipe the lenses gently.

Protect the device from dust and moisture. Store it in the supplied bag or transportation packaging. The batteries should be removed from the unit if it has not been used for a long time.

Troubleshooting

Error	Solution
No picture visible (when observing with your eyes)	<ul style="list-style-type: none"> • switch light on • adjust focus
Camera picture flickers	<ul style="list-style-type: none"> • if necessary, adjust resolution of the video graphics board
Software Installation reports „not XP approved“	<ul style="list-style-type: none"> • confirm by clicking <OK>

Specifications

System requirements for microscope camera

Minimum system requirements: PC with an Intel Pentium IV processor or higher; Windows XP with Service Pack 3*, Windows Vista (32/64Bit) with Service Pack 2* or Windows 7 (32/64Bit) with Service Pack 1*; .NET Framework 4.0*; min. 1024 MB RAM (64Bit = 2048 MB); min. 500 MB free hard drive space; free USB port; CD/DVD/BD drive.

*available for free via Windows Update (Internet connection needed)

Magnification table

Eye p.	Object.	Magn.	with Barlow I.
5x	4x	20x	40x
5x	10x	50x	100x
5x	40x	200x	400x

Disposal

Dispose of the packaging materials properly, according to their type (paper, cardboard, etc). Contact your local waste disposal service or environmental authority for information on the proper disposal.


Do not dispose of electronic devices in the household waste!

As per the Directive 2002/96/EC of the European Parliament on waste electrical and electronic equipment and its adaptation into German law, used electronic devices must be collected separately and recycled in an environmentally friendly manner.

In accordance with the regulations concerning batteries and rechargeable batteries, disposing of them in the normal household waste is explicitly forbidden. Please make sure to dispose of your used batteries as required by law – at a local collection point or in the retail market. Disposal in domestic waste violates the Battery Directive.

Batteries that contain toxins are marked with a sign and a chemical symbol. „Cd“ = cadmium, „Hg“ = mercury, „Pb“ = lead.

Cd¹

Hg²

Pb³

- ¹ battery contains cadmium
- ² battery contains mercury
- ³ battery contains lead

CE Declaration of Conformity

Bresser GmbH has issued a „Declaration of Conformity“ in accordance with applicable guidelines and corresponding standards. This can be viewed any time upon request.

WARRANTY

The period of warranty is 2 years, beginning on the day of purchase. Please keep the cash receipt as evidence of purchase. Devices which become defective during the warranty period can be returned to the dealer where the device was bought. The repaired device or a new one will then be returned to you. In the case of defects which occur after the end of the warranty period, the devices can also be returned. However, repairs which become necessary after the end of the warranty period will be subject to a service fee.

Important:

Make sure to return the device carefully packed in the original packaging in order to prevent transport damage. Please also enclose the cash receipt (or a copy). This warranty does not imply any restriction of your statutory rights.

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