- DE Experimente**
- Praktische Anleitungen
- EN Experiments**
- Practical Instructions
- FR Expériences**
- Notice pratique
- NL Experimenten**
- Gebruikssuggesties
- IT Esperimenti**
- Istruzioni pratiche
- ES Experimentos**
- Instrucciones prácticas
- PT Experiências**
- Instruções práticas

Praktische Anleitungen

Einführung

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Practical Instructions

Introduction

I would like to give you a few tips about how to take a better look at the wonderful world of microorganisms and crystals. I'll tell you how to prepare your object, for example, so that you can observe it with the microscope. The numerous experiments described should make you curious and want to use your microscope more.

Experiments with the Microscope

Experiment No. 1:

Black and White Print

Objects:

1. a small piece of paper from a newspaper with a black and white picture and some text,
2. a similar piece of paper from a magazine.

Place both pieces of paper next to one another on the microscope stage

and set your microscope to the lowest magnification. Take a look through the eyepieces and compare:

The letters on the newspaper look frayed and broken, since they are printed on raw, low-quality paper. The letters on the magazine look smoother and more complete. The pictures in the newspaper are made up of many tiny dots, which appear slightly smudgy. The pixels (halftone dots) of the magazine picture are clearly defined.

Experiment No. 2: Color Print

Objects:

1. a small piece of color printed newspaper,
2. a similar piece of paper from a magazine.

Place both pieces of paper next to one another on the microscope stage and set your microscope to the lowest magnification. Take a look through the eyepieces and compare:

The colored pixels of the newspaper often overlap. Sometimes, you'll even notice two colors in one pixel. In the magazine, the dots appear clear and rich in contrast. Look at the different sizes of the pixels.

Experiment No. 3: Textile fibers

Objects and accessories:

1. threads from various fabrics (e.g. cotton, linen, sheep's wool, silk, rayon, etc.),
2. two needles.

Place the different threads on a table and use the needles to fray them a bit. Dampen the threads with a little water. Set your microscope to the lowest magnification and observe the different textile threads through the eyepieces. Compare:

Cotton fibers come from a plant, and look like a flat, twisted ribbon under the microscope. The fibers are thicker and rounder at the edges than in the middle. Cotton fibers are basically

long, collapsed tubes. Linen fibers also come from a plant, and they are round and run in one direction. The fibers shine like silk and exhibit countless bulges on the thread. Silk comes from an animal and is made up of solid fibers that are small in diameter, in contrast to the hollow plant-based fibers. Each fiber is smooth and even and looks like a tiny glass tube. The fibers of the sheep's wool also come from an animal. The surface is made of overlapping sleeves that look broken and wavy. If possible, compare sheep's wool from different weaving mills. In doing so, take a look at the different appearance of the fibers. Experts can determine which country the wool came from by doing this. Rayon is a synthetic material that is produced by a long chemical process. All the fibers have solid, dark lines on the smooth, shiny surface. After drying, the fibers curl into the same position. Observe the differences and the similarities.

Experiment No. 4: Table Salt

Object: normal table salt.

Place a piece of black paper on the microscope stage. Then, put a few grains of salt on the paper and set your microscope to the lowest magnification. Observe the salt crystals through the eyepieces of your microscope and notice:

The crystals look like tiny dice and all have the same shape.

Experiment No. 5: Leaves and Needles

Object: 3-4 different leaves or needles from deciduous trees or fir trees.

When you go for a walk in the forest with your parents, you can collect different types of leaves and needles.

At home, place them next to one another on the microscope stage. Take a look through the eyepieces and observe the different appearance of the

leaves and needles with the lowest magnification.

Observe:

The leaves of the deciduous trees have different but more or less regular sections that are separated by lines. These are called "cells." Most often, the underside of the leaf looks different than the top, and the color is brighter. The stalk of the leaf runs through the middle. At its thicker end, there is a "lump" with a bulge. That is the part that connected the leaf to the tree, before it fell away. Some leaves also have a stalk upon which multiple leaves grow from other stalks.

Fir needles are long, thin and round. Like the leaves of deciduous trees, they have a light bulge on one side, where they grew from the tree. They do not have individual "cells," however, but look like they grew in one part. However, when you look more closely, you can see that the needle has many sections. These sections come from the step-by-step growth of the needles.

Experiment No. 6: Stones and Crystals

Object:

Stones and crystals (accessories).

You microscope comes with a few stones and crystals. You can place them next to each other on the microscope stage and than observe them through the eyepieces.

You will notice that not only the colors of the individual rock samples are different, but also the shape.

Tip:

You can also create an observation diary. In it, you can record your observations, draw pictures of particularly interesting objects and make notes.

In this way, you can look at many more objects, such as small organisms (flies, spiders, etc.) or other things from your daily life. Simply place everything on the microscope stage and take a look through the eyepieces.

You can discover so many things that you did not know before. Just give it a try!

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Bresser GmbH